



PAPER POINTS



PP POINTS



PP ONE POINTS



PP SUPER POINTS

Conformity / Regulatory Information

Generic Device Group	Absorbent Paper Points				
Model	PP Points, PP ONE Points, PP SUPER Points				
Regulation	Medical Device Regulation (EU) 2017/745				
Classification and Rule	Class I, Rule 5				
Basic UDI-DI	697060203APP0016K				
Nomenclature	<table> <tr> <td>EMDN</td><td>Q01010401</td></tr> <tr> <td>GMDN</td><td>38777</td></tr> </table>	EMDN	Q01010401	GMDN	38777
EMDN	Q01010401				
GMDN	38777				
Intended Use	Absorbent Paper Point is used for drying, cleaning, and absorbing fluids in the root canal during endodontic treatment.				
Indications of use	- For professional dental application only - Non sterile - Single use				
Other	- 28 mm length - ISO colour-coded sizes for easy identification - High absorbency for efficient canal drying - Pre-measured depth markings (16, 18, 19, 20, 22, 24 mm)				

Product Specification

Size and Specification	Model	Size	Taper (%)	
	PP Points	15 20 25 30 35 40 45 50 55 60 70 80 Assorted 15-40 or 45-80	2	
		15 20 25 30 35 40 45 Assorted 15-40 or 45-80	4	
		15 20 25 30 35 40 45 Assorted 15-40 or 45-80	6	
	PP ONE Points	R25	6	
		R40	4	
		R50	4	
		Assorted R25-R50	-	
	PP SUPER Points	F1	6	
		F2	6	
		F3	6	
		F4	5	
		F5	4	
		Assorted F1-F3	-	

Packaging Description	PP Points	Taper 2%	200 pcs/pack
		Taper 4%	100 pcs/pack
		Taper 6%	100 pcs/pack
	PP ONE Points PP SUPER Points		100 pcs/pack
Composition	Pulp, cellulose		
Traceability	UDI / LOT		
Storage Conditions	- Do not store in direct sunlight, store in a cool dry place - Temperature: 5 – 30°C		

Manufacturer Information

Manufacturer	Chengdu Sani Medical Equipment Co., Ltd.	Manufacturer SRN	CN-MF-000011715
EU Representative	Riomavix S.L.	EU Representative SRN	ES-AR-000001202
EU Importer	Premium Plus Poland sp. z o.o.	Country of Manufacture	China