

2202 patient bibs

3-ply embossed patient bibs

Two absorbent paper layers and one waterproof backing layer for maximum protection and patient comfort.


Reliable barrier against fluids and splatter, designed for daily clinical use.

product specifications


- ◆ Size: 45 × 33 cm
- ◆ 3-ply embossed construction
- ◆ Convenient pack size: 50 pcs/pack
- ◆ Single use



3-layer protection



waterproof PE
backing to protect
patient clothing



2 absorbent paper
layers for effective
liquid control




available in 14 colours


A wide range of colours to suit every practice style, branding or preference.




dark blue



blue




jade green



light green




fairy teeth pattern



lime



yellow




orange



pink




purple



burgundy



black



white



animals pattern

packaging

MOQ: 1 carton (10packs)



carton:
10 x 50 pcs

pack:
50 pcs

Order Code	Product Description	Carton Quantity	Carton Weight	Cartons per EU Palette	Carton Dimensions
2202	Patient Bibs	10 x 50 pcs/pack	3.5 kg	64 cartons	36 x 30 x 22 cm